GUIDE™ XT

Software User Guide

Version 2.0.2

Directions for Use	en
Instrucciones de uso	es
Mode d'emploi	fr
Gebrauchsanweisung	de
Istruzioni per l'uso	it
Gebruiksaanwijzing	nl
Bruksanvisning	sv
Käyttöohje	fi
Bruksanvisning	no
Brugsanvisning	da
Instruções de Utilização	pt
Instruções de uso	ptBR
Kullanım Talimatları	tr
Указания по использованию	ru
إرشادات الاستخدام	ar
Οδηγίες χρήσης	el
Návod k použití	cz
Návod na použitie	sk
Instrukcja obsługi	pl
Használati utasítások	hu

en

This page intentionally left blank

TABLE OF CONTENTS

GENERAL INFORMATION 1

Contact Data and Legal Information 1

 Contact Data 1

 Legal Information 2

 Version 2

Symbols 3

 Symbols Used in This Guide 3

Intended Use 4

 Using the System 4

Compatibility with Medical Devices 5

 Medical Software 5

 Non-Brainlab/Boston Scientific Software 5

Technical Specifications 6

 Hardware Requirements 6

 Operating Systems 6

Documentation 7

 Overview 7

 Used Abbreviations 7

 Additional Information 7

USING GUIDE XT 8

Starting a Session 8

Screen Layout 10

 Layouts and Views 10

 Toolbar 14

 Simulation Toolbar 15

 Data Menu 16

Saving Data and Closing the Software 22

INDEX 24

Troubleshooting 25

en

This page intentionally left blank

1 GENERAL INFORMATION

1.1 Contact Data and Legal Information

1.1.1 Contact Data

Legal Manufacturer
Boston Scientific Neuromodulation Corporation
25155 Rye Canyon Loop
Valencia, CA 91355 USA
(866) 789-5899 in US and Canada
(661) 949-4000, (661) 949-4022 Fax
(866) 789-6364 TTY
www.bostonscientific.com
Email: neuro.info@bsci.com

1.1.2 Legal Information

Trademarks

Microsoft® and Windows® are registered trademarks of Microsoft Corporation. All trademarks are the property of their respective holders.

Integrated 3rd-Party Software

- This software is based in part on Xerces C++ 3.1.1, developed by the Apache Software Foundation. For a full description of copyrights, disclaimers and license see <http://xerces.apache.org/>.
- This product includes libtiff 4.0.4beta. For a full description of copyrights, disclaimers and license see <http://www.remotesensing.org/libtiff/>.
- This product includes software developed by Telerik, Inc., Copyright © 2002-2016
- This software is based in part on the work of the Independent JPEG Group.
- This software contains the OpenJPEG library, which was imported and enhanced by libjpeg-turbo. For a full description of copyrights, disclaimers and license see <http://www.openjpeg.org> and <http://www.libjpeg-turbo.org>.
- Portions of this software are based on the work of Sun Microsystems Inc.

1.1.3 Version

Please verify that the installed software matches the version in this manual - GUIDE XT 2.0.2. Use the DB-7106-G2 GUIDE XT Installer to install this software on the Boston Scientific Clinician Programmer (DB-7161 or NM-7161).

1.2 Symbols

1.2.1 Symbols Used in This Guide

Warnings

Warnings are indicated by triangular warning symbols. They contain safety-critical information regarding possible injury, death or other serious consequences associated with equipment misuse.

Cautions

Cautions are indicated by circular caution symbols. They contain safety-critical information regarding possible problems with the device. Such problems include device malfunctions, device failure, damage to device or damage to property.

Notes

NOTE: Notes are formatted in italic type and indicate additional useful hints.

1.3 Intended Use

1.3.1 Using the System

Guide XT Intended Use

Guide XT is intended to be used post-operatively in an office environment on a computer as a planning tool for the programming of the Boston Scientific DBS system.

Indications for Use

Guide XT is intended to display medical images and simulate stimulation output. It includes functions for image manipulation and 3D visualization (reconstructions and volume rendering). Features include the display of a simulated DBS lead from a patient's CT scan compared to an anatomical atlas.

Intended User

The intended users are healthcare professionals, educated for the planning and execution of DBS procedures. These are, in general, neurosurgeons and neurologists. Usually the overall procedure is split between surgical and neurological responsibilities and are performed accordingly.

Place of Use

The intended use environment is a (laptop) computer in a hospital (in office and operating room environment) or a doctor's office. The common use environment will be an office environment.

Careful Handling

Only trained medical professionals may operate system components and accessory instrumentation.

Plausibility Review

Before patient treatment, review the plausibility of all information input to and output from the system.

Responsibility

This system solely provides additional assistance to the medical professional and does not by any means substitute or replace the medical professional's experience and/or responsibility during its use.

1.4 Compatibility with Medical Devices

1.4.1 Medical Software

Compatible Medical Software

Guide XT is compatible with:

- Brainlab **Content Manager** 2.3 and 2.4
- Brainlab **Patient Selection** 4.1
- Brainlab **DICOM Viewer** 2.2, 3.1, and 3.2
- Brainlab **Image Fusion** 2.0 and 3.0
- Brainlab **Smartbrush** 2.1, 2.5, and 3.0
- Brainlab **Fibertracking** 1.0
- Brainlab **Anatomical Mapping** 1.0
- Brainlab **Object Manipulation** 1.0

Other Brainlab/Boston Scientific Software

Other compatible Brainlab/Boston Scientific software may become available after the release of this user guide. If you have questions regarding compatibility of software contact Boston Scientific support. If you are running software versions other than those specified above, contact Boston Scientific support for clarification regarding compatibility with those devices.

1.4.2 Non-Brainlab/Boston Scientific Software

Authorization

Only authorized Brainlab or Boston Scientific employees may install software on the Brainlab system. Do not install or remove any software applications.

Compatible Non-Brainlab Software

Guide XT is compatible with:

- Microsoft Surface Pro 3
- Microsoft Windows 7
- Microsoft Windows 8

For detailed and up-to-date information regarding compatible operating systems, please contact Boston Scientific support.

Other Non-Brainlab/Boston Scientific Software

Only software specified by Brainlab or Boston Scientific may be installed and used with Guide XT.

1.5 Technical Specifications

1.5.1 Hardware Requirements

	Minimum
Processor	4 logical cores (e.g. Intel Core i5 with or comparable processor)
Memory	4 GB
Graphics Card	DirectX 10.1 compatible 512 MB graphics memory
Screen Resolution	1280 x 1024
Color Depth	24bit
Disk Space	20 GB

1.5.2 Operating Systems

The elements are released only for the following desktop systems (server systems not listed here):

- win7-x64-sp1
- win8.1-x64
- win10-x64

1.6 Documentation

1.6.1 Overview

Intended Audience

This user guide is intended for medical professionals and/or their staff.

Reading User Guides

The user guides describe complex medical devices and software that must be used with care.

It is important that all users of system, instruments and software:

- Read the user guides carefully before handling the equipment
- Have access to the user guides at all times

1.6.2 Used Abbreviations

Abbreviations

This user guide may contain the following abbreviations:

Abbreviation	Definition
DBS	Deep brain stimulation
IPG	Implantable pulse generator
SFM	Stimulation field model

1.6.3 Additional Information

For indications, contraindications, warnings, precautions, adverse events summary, sterilization, component disposal, storage and handling, and warranty information, refer to the Information for Prescribers DFU for your Vercise™ DBS System. For other device specific information not included in this manual or labeling symbols, refer to the appropriate DFU for your Vercise™ DBS System as listed on your DBS Reference Guide.

GUIDE XT supports all Boston Scientific, CE marked, DBS Implantable Pulse Generators.

2 USING GUIDE XT

2.1 Starting a Session

General Information

Guide XT is designed to support neurosurgeons and neurologists in Deep Brain Stimulation (DBS) treatments. It includes stimulation field simulation to aid in appropriate DBS parameter adjustment. It does not directly interact with the DBS programming equipment.

Typical Workflow

The typical **Guide XT** workflow includes the use of both required and optional additional applications. For more information, see the **Software User Guides** for these applications.


Figure 1

Steps
1. Start the Guide XT workflow.
2. Select patient data in Patient Selection .
3. Optional: Fuse image sets using Image Fusion .
4. Optional: Segment structures using Anatomical Mapping .
5. Localize and define leads using Lead Localization .
6. Select the Guide XT button in Content Manager .

Compatible Lead Types

The following Boston Scientific lead types are compatible with **Guide XT**:

- DB2201-Standard
- DB2202-Directional

Configuration Dialog

When **Guide XT** is started, the configuration dialog opens. Here you must define the IPG used and assign the leads that were defined in **Lead Localization** to the corresponding IPG contacts.

*NOTE: The configuration dialog can also be accessed later via the **Configuration** button on the main screen.*


Figure 2

No.	Component
①	Model: Select IPG type.
②	Implant Date: Select implantation date.
③	Battery Type and Model Number: Displayed according to IPG type.
④	Patient information
⑤	Assign lead(s) to the IPG contacts to which it is connected. <i>NOTE: Only valid IPG-lead configurations can be assigned.</i>

2.2 Screen Layout

2.2.1 Layouts and Views

Standard Layout

The **Standard** layout combines a large view that shows only 3D structures with an **Inline** and a **Perpendicular** view that show additional image information. Only basic rotation and the **Zoom** functions are available. View manipulation is performed in-view.



Figure 3

No.	Component
①	Views: <ul style="list-style-type: none"> • 3D view • Inline view: The yellow depth line is adjustable. Changes are reflected in the Perpendicular view. • Perpendicular view
②	Simulation toolbar
③	Toolbar

Overview Layout

The **Overview** layout combines a large view with full visualization and interaction functionality with an **Inline** and a **Perpendicular** view.

No.	Component
①	<p>Views:</p> <ul style="list-style-type: none"> • 3D view: <ul style="list-style-type: none"> - Switch between Axial, Coronal and Sagittal views. - Select Aligned to align the view along the trajectory of the lead. • Inline view: The yellow depth line is adjustable. Changes are reflected in the Perpendicular view. • Perpendicular view <p><i>NOTE: Use the Rotate function to dynamically switch between Axial, Coronal and Sagittal views.</i></p>
②	Simulation toolbar
③	Toolbar


How to Switch Screen Layouts

Figure 4

Step
In the Data menu, select Standard or Overview layout.

View Interaction

In the **Inline** and **Perpendicular** views:

- Rotation made in one view, is reflected in the other.
- Objects are shown in 2D.
- The medial/lateral directions are identified in the **Perpendicular** view.



In the **3D** view:

- The dot ① represents the intersection of the yellow line ② with the lead in the **Inline** view
- The ruler ④ is only visible in the **3D** view when you are in-plane with axial, coronal, sagittal orientations.

For directional leads, the orientation is represented by a directional marker ③.

2.2.2 Toolbar

Main Menu Functions

Button	Function
 Scroll	Scrolls in 0.5 mm increments.
 Zoom	Zooms to the region of interest.
 Pan	Pans vertical and horizontal planes in the 3D view.
 Rotate	Rotates the view.
 Windowing	Slide left or right across the viewing area to adjust the contrast and slide up or down to adjust the brightness.
 Center	Centers the view to the center of the electrode contacts.
 Configuration	Opens the Configuration dialog (see page 9).
 Create Report	Creates treatment report (see page 22).

Active Lead

When a lead is selected the first time, all objects assigned to the opposite hemisphere are initially hidden. For example, when selecting a lead located in the left hemisphere, all objects assigned to the right hemisphere are hidden in the views.

How to Select Active Lead

Figure 5

Step
Toggle to select active lead ①.

2.2.3 Simulation Toolbar

Simulation Parameters

Figure 6

No.	Component
①	Simulation parameter label
②	Create button: Select to create a new simulation.
③	Copy button: Select to create a new simulation parameter that copies the current simulation parameters.

Comment Field and Colors

Figure 7

No.	Component
①	Select the Comment field to enter an optional comment.
②	Select a color to assign it to the simulation parameter label.

2.2.4 Data Menu

How to Access Data Menu

Step	
	Select Data button. The Data menu pops out.

Selected Plan


Figure 8

The selected plan status and name ① is displayed within the **Data** menu.

Data Menu: Images


Figure 9

In the **Images** drop down, you can select the images to display in the views.

Data Menu: Objects

Objects may include auto-segmented or drawn objects.


Figure 10

No.	Component
①	Toggle 2D/3D objects (in 3D view)
②	Show/hide all objects (in 3D view) <i>NOTE: You can also show/hide individual objects. When outlined in blue, an object is set to visible.</i>
③	Object characteristics: <ul style="list-style-type: none"> • Name assigned to the object • Type • Comment
④	Color selection
⑤	Delete object

Data Menu: Fiber Bundles


Figure 11

No.	Component
①	Toggle 2D/3D fiber bundles (in 3D view)
②	Show/hide fiber bundles

Data Menu: Leads


Figure 12

No.	Component
①	Show/hide Leads
②	Lead characteristics: <ul style="list-style-type: none"> • Lead Type • Hemisphere that is assigned to the lead • Implant Date • Comment

Data Menu: Simulation Parameters


Figure 13

No.	Component
①	<ul style="list-style-type: none"> • Basis: The lead to which the parameters are assigned. • Name: assigned to the parameters. • Comment field: Review and/or edit comments entered on the main screen.
②	Color: Assign a color identifier.
③	Export: Export as DICOM object. The exported object is then visible in the Objects dropdown and can be used in other applications that support DICOM objects.
④	Delete: Delete the parameters.

NOTE: Only parameters for the currently selected lead are shown. The selected parameters are highlighted orange.

en

2.3 Saving Data and Closing the Software

How to Create a Report


Figure 14

Steps	
1.	 Select Create Report .
2.	In the Save Plan dialog, edit the Plan Name if needed.
3.	<ul style="list-style-type: none"> • Select Save to save current patient data and proceed to report. • Select Cancel to proceed to report without saving plan.
4.	Report is created, including: <ul style="list-style-type: none"> • Patient Information • Hardware Information • IPG Simulation Parameters • Screenshots <p>In the preview dialog, you can page through, save or print the report.</p>

Saving and Closing Guide XT

Button	Explanation
 Back	Returns you to the previous step and discard changes.
 Done	Select Done to perform the following: <ul style="list-style-type: none">• Save your existing task• Close Guide XT• Open another application
 Home	Pressing Home returns you to Content Manager and saves changes.

INDEX

	A			O
active lead		18	objects	21
	B		overview layout.....	14
back.....		26		P
	C		pan.....	17
center.....		17	PDF report.....	25
change layout		15		R
close software.....		26	rotate	17
compatible lead		11		S
configuration.....		12,17	scroll	17
create report.....		17,25	select lead	18
	D		simulation parameter.....	18,24
data menu.....		19	standard layout.....	13
documentation.....		10	starting a session.....	11
done.....		26	switch layout.....	15
	F			T
fiber bundles.....		22	treatment report.....	25
	H			V
home.....		26	view layout.....	15
	I		views.....	16
images.....		20		W
image selection.....		20	windowing.....	17
image view.....		16		Z
	L		zoom.....	17
leads.....		23		
lead selection.....		23		
	M			
menu buttons.....		17		

Troubleshooting

There are no user serviceable parts. If you have a specific question or issue and need to contact Boston Scientific, choose your locality from the following list:

Argentina

T: +5411 4896 8556 F: +5411 4896 8550

Australia / New Zealand

T: 1800 676 133 F: 1800 836 666

Austria

T: +43 1 60 810 F: +43 1 60 810 60

Balkans

T: 0030 210 95 37 890 F: 0030 210 95 79 836

Belgium

T: 080094 494 F: 080093 343

Brazil

T: +55 11 5853 2244 F: +55 11 5853 2663

Bulgaria

T: +359 2 986 50 48 F: +359 2 986 57 09

Canada

T: +1 888 359 9691 F: +1 888 575 7396

Chile

T: +562 445 4904 F: +562 445 4915

China – Beijing

T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou

T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai

T: +86 21 6391 5600 F: +86 21 6391 5100

Colombia

T: +57 1 629 5045 F: +57 1 629 5082

Czech Republic

T: +420 2 3536 2911 F: +420 2 3536 4334

Denmark

T: 80 30 80 02 F: 80 30 80 05

Finland

T: 020 762 88 82 F: 020 762 88 83

France

T: +33(0)1 39 30 97 00 F: +33(0)1 39 30 97 99

Germany

T: 0800 072 3301 F: 0800 072 3319

Greece

T: +30 210 95 42401 F: +30 210 95 42420

Hong Kong

T: +852 2960 7100 F: +852 2563 5276

Hungary

T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore

T: +91 80 5112 1104/5 F: +91 80 5112 1106

India – Chennai

T: +91 44 2648 0318 F: +91 44 2641 4695

India – Delhi

T: +91 11 2618 0445/6 F: +91 11 2618 1024

India – Mumbai

T: +91 22 5677 8844 F: +91 22 2617 2783

Italy

T: +39 010 60 60 1 F: +39 010 60 60 200

Korea

T: +82 2 3476 2121 F: +82 2 3476 1776

Malaysia

T: +60 3 7957 4266 F: +60 3 7957 4866

Mexico

T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East / Gulf / North Africa

T: +961 1 805 282 F: +961 1 805 445

The Netherlands

T: +31 30 602 5555 F: +31 30 602 5560

Norway

T: 800 104 04 F: 800 101 90

Philippines

T: +63 2 687 3239 F: +63 2 687 3047

Poland

T: +48 22 435 1414 F: +48 22 435 1410

Portugal

T: +351 21 3801243 F: +351 21 3801240

Singapore

T: +65 6418 8888 F: +65 6418 8899

South Africa

T: +27 11 840 8600 F: +27 11 463 6077

Spain

T: +34 901 11 12 15 F: +34 902 26 78 66

Sweden

T: 020 65 25 30 F: 020 55 25 35

Switzerland

T: 0800 826 786 F: 0800 826 787

Taiwan

T: +886 2 2747 7278 F: +886 2 2747 7270

Thailand

T: +66 2 2654 3810 F: +66 2 2654 3818

Turkey – Istanbul

T: +90 216 464 3666 F: +90 216 464 3677

Uruguay

T: +59 82 900 6212 F: +59 82 900 6212

UK & Eire

T: +44 844 800 4512 F: +44 844 800 4513

Venezuela

T: +58 212 959 8106 F: +58 212 959 5328

Note: *Phone numbers and fax numbers may change. For the most current contact information, please refer to our website at <http://www.bostonscientific-international.com/> or write to the following address:*

Boston Scientific Neuromodulation Corporation
25155 Rye Canyon Loop
Valencia, CA 91355 USA

This page intentionally left blank

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- Este producto contiene software desarrollado por Telerik, Inc., Copyright © 2002-2016
- Ce produit comprend un logiciel développé par Telerik, Inc., Copyright © 2002-2016.
- Dieses Produkt enthält Software, die von Telerik, Inc. entwickelt wurde, Copyright © 2002-2016
- Questo prodotto include software sviluppato da Telerik, Inc., Copyright © 2002-2016
- Dit product bevat software ontwikkeld door Telerik, Inc., Copyright © 2002-2016
- Denna produkt innefattar programvara som utvecklats av Telerik, Inc., Copyright © 2002–2016