# Medtronic

DBS

37086

Extension Kit for Deep Brain Stimulation (8-4)

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! USA Rx only

**CE 0123**  
2011

## Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

Open here

Extension length

Do not reuse

Do not re-sterilize

Do not use if package is damaged

Sterilized using ethylene oxide

Caution, consult accompanying documents

Consult instructions for use

Date of manufacture

Manufacturer

Use by

Serial number

**CE 0123**  
2011

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

Authorized representative in the European community

For USA audiences only

Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

# Table of contents

**Description 7**

**Package contents 7**

**Device specifications 7**

**Instructions for use 10**

Preparing the neurostimulator implantation site 10

Creating a subcutaneous pocket for the neurostimulator 10

Removing the lead cap 10

Tunneling the extension 13

Connecting the lead to the extension 14

Removing the percutaneous extension (if applicable) 16

## **Information available for the Deep Brain Stimulation Therapies:**

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, and component disposal. The information for prescribers addendums contain indication-specific information. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility and battery longevity manual describes programming considerations and provides battery longevity information to aid in the appropriate neurostimulator selection.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

**USA** The clinical summary provides information about the clinical study results for the neurostimulation system.

## Description


The Medtronic Model 37086 Extension is part of a neurostimulation system for deep brain stimulation.

## Package contents

- Extension
- Connector boots (4)
  - cylindrical (2)
  - winged (2)
- Octapolar in-line neurostimulator plug
- Extra setscrews (2)
- Torque wrench
- Product literature
- Warranty card (USA only)
- Registration form (USA only)

## Device specifications

The extension has four connectors on the distal end and eight connectors on the proximal end. The in-line distal end connects to a Medtronic lead. The proximal end connects to a Medtronic neurostimulator.

**Figure 1.** Model 37086 extension.

**Table 1. Device specifications for the Model 37086 Extension<sup>a</sup>**

<b>Description</b>	<b>Model 37086</b>
Resistance <sup>b,c</sup>	Maximum 38.0 $\Omega$ for all lengths
Length	40 cm, 60 cm, or 95 cm
Distal end (lead)	
Connector	Quadripolar, in-line
Contact spacing	4.3 mm
Lead entrance diameter	1.5 mm
Outer diameter	3.8 mm
Proximal end (neurostimulator)	
Connector	Octapolar, in-line
Contact spacing	2.8 mm
Diameter	1.3 mm

<sup>a</sup> All measurements are approximate.

<sup>b</sup> Electrical resistance of this device only.

<sup>c</sup> Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

**Table 2. Material of components in the Model 37086 package**

<b>Component</b>	<b>Material</b>	<b>Material contacts human tissue</b>
Extension		
Conductor wire	Silver core MP35N	No
Distal end (lead)		
Insulation	Silicone rubber	Yes
Setscrew connector block	Stainless steel	Yes
Proximal end (neurostimulator)		
Conductor wire insulation	Fluoropolymer	No
Contacts	MP35N	Yes
Insulation	Polyurethane	Yes
Connector boots		
Clear	Silicone rubber	Yes
White radiopaque	Silicone rubber, barium sulfate	Yes
Neurostimulator plug	Polyurethane	Yes
Setscrews	Titanium	Yes

## Instructions for use

Implanting physicians should be experienced in stereotactic and functional neurosurgery and deep brain stimulation procedures, as well as thoroughly familiar with all product labeling.

### **Cautions:**

- Do not bend, kink, stretch, or twist the lead or extension, which may damage the component. Component damage would require surgical replacement.
- Do not use any instrument to handle the extension. The force may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.

## Preparing the neurostimulator implantation site

1. Select a location for the neurostimulator and position the patient accordingly.

**Caution:** Select a location that meets the following criteria:

- A minimum of 20 cm (8 in) away from another Medtronic neurostimulator to minimize telemetry interference and possible inappropriate therapy
- On the opposite side of the body from another active implanted device (for example, a pacemaker or defibrillator) to minimize possible interaction between the devices
- Away from bony structures to minimize discomfort at the neurostimulator site
- Away from areas of restriction or pressure to minimize the potential for skin erosion and patient discomfort
- Accessible to the patient for proper neurostimulator recharging and proper operation of a patient control device

**Note:** Also take into account:

- Possible future cardiac needs (eg, pacemaker, defibrillator). Implanting a neurostimulator on the patient's right side allows for possible future placement of cardiac devices on the patient's left side.
- The cosmetic needs of the patient.

2. Identify the neurostimulator pocket site and mark the intended extension route from the lead-incision site to the neurostimulator pocket.

## Creating a subcutaneous pocket for the neurostimulator

1. Make an incision equal to the length of the neurostimulator.
2. Use blunt dissection to create a subcutaneous pocket.

**Note:** Refer to the neurostimulator implant manual for the proper subcutaneous pocket depth.

## Removing the lead cap

**Note:** Before opening the extension package, verify the model number, use-by date, extension length, and connector type.

For each implanted lead:


1. Locate the lead cap at the proximal end of the lead.

**Note:** Use fluoroscopic observation or palpate the area with the fingers to locate the lead cap and confirm access without risk of cutting the lead.

2. Make an incision sufficient to expose the lead and lead cap. Allow room to hold the lead firmly when accessing the lead cap and when connecting the lead to the extension (Figure 2).

**Cautions:**

- Be extremely careful when using sharp instruments around the lead body to avoid nicking or damaging the lead, which would require surgical replacement.
- Do not pull the lead cap out of the tissue to prevent excess stress on the lead. If more slack is needed, gently grasp the lead body and lift the lead cap out of the tissue.
- Do not bend, kink, stretch, or twist the lead body when accessing the lead cap. Lead damage may occur, which would require surgical replacement.


**Figure 2.** Make incision to expose lead and lead cap.

3. Cut the distal ligature over the lead end of the connector boot (Figure 3).


**Caution:** Do not cut the lead when removing the ligature. Cutting the lead insulation can result in loss of stimulation and lead failure, which would require surgical replacement.


**Figure 3.** Cut distal ligature over lead end of connector boot.

4. Cut the connector boot over the lead cap to expose the setscrew (Figure 4).

**Figure 4.** Cut connector boot to expose setscrew.

5. Hold the setscrew connector block firmly between thumb and two fingers to prevent rotation (Figure 5).

**Figure 5.** Hold setscrew connector block between thumb and two fingers.


- Using the torque wrench, loosen the setscrew in the setscrew connector block by turning the wrench counterclockwise (approximately one turn) (Figure 6).

**Note:** Ensure the torque wrench is perpendicular to the setscrew surface while loosening the setscrew.

**Cautions:**

- Ensure the torque wrench is perpendicular to the setscrew surface while loosening the setscrew. Failure to hold the wrench perpendicular while loosening the setscrew may result in damage to the lead, resulting in intermittent or loss of stimulation.
- Hold the setscrew connector block firmly between thumb and two fingers to support the lead body and prevent rotation. Failure to properly prevent rotation may break or damage the lead, and may require lead replacement.


**Figure 6.** Loosen setscrew.

- Stabilize the lead and gently remove the lead from the setscrew connector block (Figure 7).


**Caution:** If resistance is felt while removing the lead from the lead cap, loosen the setscrew slightly to ensure that it clears the lead contact. Avoid disengaging the setscrew. Inspect the lead for damage (flattening or stretching of the lead) if resistance was felt prior to removal, which would require surgical replacement.

**Figure 7.** Remove lead from setscrew connector block.

8. Withdraw the lead cap and connector boot through the incision and discard (Figure 8).

**Figure 8.** Withdraw lead cap and connector boot.


## Tunneling the extension

**Warning:** Do not place the lead-extension connector in the soft tissues of the neck. Placement in this location has been associated with an increased incidence of lead fracture, which would require surgical replacement.

**Warning:** Avoid tunneling too deeply or superficially. Injury to nerve tissue (such as the spinal accessory nerve), vascular injury, and tunneling through unintended anatomy can occur, resulting in serious injury, prolonged hospitalization, or surgical intervention.

**Cautions:**

- When routing the extension, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When two leads are implanted, route the lead-extensions so the area between them is minimized (Figure 9). If the lead-extensions are routed in a loop, the loop will increase the potential for electromagnetic interference (EMI).

**Figure 9.** Routing for multiple lead-extensions.

**Note:** If the procedure is not done under general anesthesia, administer the appropriate sedation to reduce patient discomfort during the tunneling portion of this procedure.

1. Identify the tunneling path. Tunnel from the lead site to an intermediate incision on either the top portion of the parietal bone or the mastoid area, to the neurostimulator pocket.
2. Make a secondary pocket at the intermediate incision site to contain the lead-extension so that the lead-extension connection does not lie directly below the incision.
3. Refer to the tunneling tool manual for detailed instructions on how to assemble the tunneling tool and perform the tunneling procedure.

## Connecting the lead to the extension

### Cautions:


- Do not use saline or other ionic fluids at connections, which could result in a short circuit.
- Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

Two different boot styles are provided, winged and cylindrical.

For each lead:


1. Slide the appropriate boot(s) over the proximal end of the lead(s):
  - Push the clear boot over the proximal end of the lead in a single-lead system or over the first lead (left hemisphere) in a dual-lead system (Figure 10A).
  - Push the white radiopaque boot over the proximal end of the second lead (right hemisphere) in a dual-lead system (Figure 10B).

**Note:** For dual-lead systems, the clear boot will distinguish the first lead-extension (left hemisphere) from the white radiopaque boot over the second lead-extension (right hemisphere).

**Figure 10.** Push the clear boot over the first lead and the white radiopaque boot over the second lead.


2. Wipe the lead body and the extension setscrew connector junction with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution to rinse the connections. After rinsing, dry all connections.
3. Hold the extension securely in your fingers. Insert the lead fully into the extension setscrew connector (Figure 11).

**Figure 11.** Insert the lead fully into the extension setscrew connector.

**Note:** Each lead contact must be aligned under each setscrew.

4. Use the torque wrench to tighten each setscrew (turning clockwise) to complete the electrical circuit. Hold the extension firmly on the sides of each setscrew connector block while tightening the setscrew. Tighten until the torque wrench clicks at least once (Figure 12).


**Figure 12.** Tighten the setscrews holding the extension firmly on the sides of each setscrew connector block.

**Cautions:**

- To avoid overtightening, use only the torque wrench included in the extension package to tighten the extension setscrews. Overtightening the extension setscrews may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.
- Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.
- Do not pull the lead body taut when implanted. The extension is available in different lengths. Select an extension length that allows connection without tension. Excess tension on the lead or extension may cause lead migration or lead or extension fracture, resulting in a loss of therapy and requiring surgical intervention.
- Ensure the torque wrench is perpendicular to the setscrew surface while tightening the setscrew. Failure to hold the wrench perpendicular while tightening the setscrew may result in undertightening or false tightening of the setscrew, resulting in intermittent or loss of stimulation.
- Hold the setscrew connector block firmly between thumb and two fingers to support the lead body and prevent rotation. Failure to properly prevent rotation may break or damage the lead or extension, and may require lead or extension replacement.

5. Slide the connector boot over the lead-extension connection (Figure 13A) and tie a ligature (using nonabsorbable silk) around both ligature channels of the boot (Figure 13B).

**Figure 13.** Slide boot over lead-extension connection and tie ligatures around the ligature channels on both ends of the boot.

**Note:** Sterile water may be used as a lubricant for ease of placing the boot.

**Cautions:**

- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure and surgical replacement.
- Do not tie ligatures around the lead or extension body, which may damage the insulation. Damaging the insulation can result in loss of stimulation and component failure, which would require surgical replacement.
- Do not overtighten ligatures on the connector boot. Ligatures that are too tight may damage the component, which would require surgical replacement.
- Ensure that the end of the connector boot is secure to prevent fluid in the lead-extension connection, which may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

- a. If using a winged boot, make sure the flat side of the boot lies against the skull.
- b. If using the cylindrical boot, suture it to the underlying tissue to minimize the possibility of migration. Place this suture in the suture channel at the proximal end of the boot covering the lead-extension connection. Loop this suture through the underlying tissue (periosteum) to stabilize the connection.

6. For a second lead and extension, repeat steps 1 through 5.

**Note:** If two extensions were implanted but only one lead will be connected, cover the unused extension with a connector boot provided in the extension product packaging or accessory kit. Tie a ligature (using nonabsorbable suture) around both ligature channels of the boot.

7. For instructions on making the extension-neurostimulator connection, refer to the appropriate neurostimulator manual.

## Removing the percutaneous extension (if applicable)


If a percutaneous extension was used for test stimulation at lead implant, follow these steps to remove the percutaneous extension.

For each implanted lead:

1. Withdraw the external segment of the percutaneous extension approximately 1 cm (0.4 in) from where it exits the skin (Figure 14).


△ **Caution:** Use gentle traction on percutaneous extension to avoid dislodging the lead.


**Figure 14.** Withdraw the external segment of the percutaneous extension.

2. Cut off and discard this section of extension.
3. Locate the percutaneous extension connector at the proximal end of the lead and make an incision to expose it. Allow room to hold the lead connection firmly to prevent dislodgement.
4. Cut the distal ligature on the lead end of the connection. Then cut the connector boot over the setscrews to expose the setscrews (Figure 15).

**Figure 15.** Cut the ligature and connector boot to expose setscrews.

5. Hold each setscrew connector block firmly between thumb and two fingers to prevent rotation while loosening setscrews (Figure 16).


**Figure 16.** Hold setscrew connector block between thumb and two fingers.

6. Using the torque wrench, loosen each of the four setscrews by turning the wrench counter-clockwise (approximately one turn) (Figure 17).

**Note:** Ensure the torque wrench is perpendicular to the setscrew surface while loosening the setscrew.


△ **Caution:** Hold the setscrew connector block firmly between thumb and two fingers to support the lead body and prevent rotation. Failure to properly prevent rotation may break or damage the lead, and may require lead replacement.

**Figure 17.** Loosen setscrews.


7. Stabilize the lead and gently remove the lead from the setscrew connector blocks (Figure 18).

 **Caution:** If resistance is felt while removing the lead from the percutaneous extension, loosen the setscrews slightly to ensure that they clear the lead contacts. Avoid disengaging the setscrews. Inspect the lead for damage (flattening or stretching of the lead) if resistance was felt prior to removal, which would require surgical replacement.

**Figure 18.** Remove the lead from setscrew connector block.

8. Withdraw the percutaneous extension and connector boot through the incision and discard the boot (Figure 19).

**Figure 19.** Withdraw the percutaneous extension and connector boot.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.